

LEVITAS<sup>®</sup> BIO

# Redefining Sample Preparation

The LeviCell Systems



# Accelerate Scientific Discovery With the Power of Levitation Technology

## Our Mission

LevitasBio® is dedicated to advancing science and human health by providing researchers with new and powerful methods for sample preparation and characterization. The LeviCell® platform is a first-in-class cell enrichment solution powered by a proprietary, label-free, and gentle Levitation Technology. Our unique solution generates high-quality samples while preserving the cells' native states, enabling researchers across multiple scientific disciplines to uncover true biology.

Increase sample viability up to 99%

Remove contaminating dead cells & debris

Run cell numbers as low as 8,000

Cut sample preparation time by 80%

## LeviCell 1.0



### Elevating Performance

The LeviCell 1.0 system provides breakthrough performance, elevating sample enrichment to new heights. Achieve viable cell enrichment with a 3-step workflow including a 20-minute walk-away run.

- User-defined stringency for purity or yield
- Gentle, label-free enrichment
- Real-time visualization and characterization of sample
- No run-to-run calibration

### LeviCell 1.0 Cartridge

- Single-use consumable
- 100% of the sample is contained in the cartridge
- Accommodates particles up to 350 µm in size
- Sterile irradiated option available



## All in a Single Run



### Viable Cell Enrichment

Label-free  
No pressure  
No cell stress



### Cell-Specific Selection

Easy to multiplex  
Ready-to-use kits  
Label-free cells of interest



### Debris Removal

Label-free  
Removed automatically  
No extra wash steps



### Sample Analysis

Visualize the full sample  
Cell characterization  
Insight to upstream steps

## LeviCell EOS

### Parallel Processing Power

The LeviCell EOS system is the next generation in high throughput sample preparation and characterization. An enhanced optical system enables the full visualization of the separation channel during the run. Characterize samples and draw conclusions during the sample preparation process.

- Flexibility to expand up to 4 EOS modules, controlled and configured independently
- Run 4 different samples or 4x the same sample
- Examine samples in real-time during the enrichment process
- Obtain consistent and robust viability enrichment


### LeviCell EOS Cartridge

- Same qualities of the LeviCell 1.0 cartridge with 4x throughput
- Full separation channel visualization
- Optically clear channels for enhanced sample characterization and visualization
- Fully sterile option, with sterile inlet well covers



# Leverage Physical Forces to Obtain Native State Biology

## The Workflow



### STEP 1: SAMPLE INPUT

Prepare the sample by adding the Levitation Agent, a paramagnetic solution. After inserting the cartridge into the system, dispense the sample into the inlet well.

### STEP 2: SAMPLE ENRICHMENT

The system automatically loads the sample into the separation channel while the sample is imaged and levitated. The system performs automated, label-free enrichment plus debris removal via levitation.

### STEP 3: SAMPLE COLLECTION



Once the cells have levitated to their final equilibrium position, a live view in the software enables the user to define how the sample will be divided (fraction of interest or waste). Separated fractions are collected into the cartridge outlet wells.

## How Levitation Technology Works

Levitation takes place due to forces exerted by the magnets on the fluid around the cells. Viable cells are impermeable to the Levitation Agent, resulting in strong levitation.

Dead and dying cells are permeable to the Levitation Agent due to their compromised membranes, so the magnetic forces are weaker as they are more similar to the surrounding medium.

- Viable cells levitate higher and form a tight uniform band
- Dead and dying cells levitate lower and form a diffuse band



## Enable Diverse Starting Samples

### Low viability samples

archived,  
frozen

### Low quality samples

FNA,  
core biopsies

### Sensitive cell types

microglia,  
neutrophils

### Large fragile samples

organoids,  
cardiomyocytes

### Complex samples

dissociated  
tumor cells

## Streamlined & Simplified Workflow for Faster Results

Combine several steps from several different platforms into one simplified and streamlined workflow with the Levitation Technology and LeviCell systems.

### CURRENT STATE




### INTEGRATED PROTOCOL




## Real-Time Visualization

The LeviCell systems' onboard camera provides the ability to see cells levitating during the run. This offers added capabilities that are unmatched by any other sample preparation technology available.


- Inspect sample quality by observing the presence of debris, cell clumps, and tissue chunks
- Use red and green fluorescence during levitation (optional)
- Observe cell population heterogeneity, behavior, or extract quantifiable metrics



Beginning equilibration



Ending equilibration



Lung dissociated tumor cells with debris

# Start With Debris-Free, Highly Viable Samples Every Time

## Gain Deeper Insights With More Viable Single Cells


Dissociated tumor cells (DTC) are precious and critical to understanding disease biology and are also often challenging samples to work with. Achieving high-quality scRNA-seq data from DTCs can be difficult. The LeviCell systems can overcome these hurdles by enriching viable cells with minimal loss and avoiding perturbation to cell profiles. Enable the best samples to acquire the best results.



Lung Cancer Dissociated Tumor Cells

### Primary Challenges

- Low viability
- Low cell number
- High debris
- Hard to obtain



Start with the highest quality DTC sample using the LeviCell

Produce efficient single-cell sequencing libraries


Obtain 3.5x more high-quality single cell data

Enrich 6x higher cancer cell representation using LeviSelect CD45+ Depletion Kit

## Viable & Cell-Specific Enrichment

The LeviSelect™ Depletion and Enrichment Kits empowers gentle enrichment of viable and cell-specific populations, which remain untouched. On the LeviCell systems, achieve dead cell and debris removal with simultaneous cell-specific enrichment in a single run.


- Integrate cell viability and cell-specific enrichment with Levitation Technology
- Enrich specific cells of interest easily and quickly
- Use custom enrichment kits for workflow flexibility



## Transform Low Quality Samples Into Sequencing Success

The LeviCell systems rescue challenging samples at low cell numbers and at low viability when other methods fail. Poor-quality samples lead to lost sequencing reads that negatively impact scRNA-seq data.


- Increase viability across multiple tissue types
- Enrich viable cells from low starting cell numbers
- Remove unaccounted-for debris



## Acquire High-Quality Data From Native-State Cells

Stressed cells are not representative of true biology. High % of mitochondrial reads per single cell indicates cells are undergoing oxidative stress. Cells enriched on the LeviCell systems have lower % of mitochondrial reads mapped.


- Preserve native-state cells
- Reduce data filtering for higher data quality, capture more relevant biological information
- Reflect true biology with more cells, genes, and transcripts





## Preserve All Cell Subtypes & Frequencies

Make efficient use of precious samples from diseased human tissue to obtain high-quality scRNA-seq results. The LeviCell systems perform viability enrichment and debris removal while preserving all cell types and cell frequencies within even the most fragile or complex samples.

- Gentle enrichment of fragile cell types
- Generate data that represents true biology
- Consistent across samples



Instrument Specifications	LeviCell 1.0 	LeviCell EOS 
Number of samples/run	1	4
Temperature control	No	Optional**
Visualization	Fixed position	Full separation channel
Fractionation analysis	From image collection	From lane 1 image collection
Cartridge	S2.3, S2.3 Sterile	EOS-4, EOS-4 Sterile
Sterile inlet covers	No	Yes
Magnetic core	Fixed	User exchangeable
Illumination	Brightfield (transmitted illumination) Ex470 nm, Em524/29 nm (eg Calcein) Ex560 nm, Em 628/33 nm (e.g. PI)	Brightfield (transmitted illumination) Ex470 nm, Em524/43 nm (eg Calcein) Ex560 nm, Em 628/65 nm (e.g. PI)
Instrument software	Experiment Manager	EOS Manager
Analysis software	Experiment Analyzer	LeviMetrics™**
Biosafety cabinet compatibility	Standard size	Large footprint or custom
Footprint (W x D x H)	47 cm x 34 cm x 24 cm 19" x 13" x 9"	44 cm x 63 cm x 46 cm 17.3" x 24.8" x 18.1"
Weight	15.3 kg (33.6 pounds)	45.4 kg (100 pounds)

\*Dependent on cell size, density  
\*\* Not yet available

Part Numbers	Product Description
<b>Instruments</b>	
1000001	LeviCell 1.0 with 1 year warranty
1000021	LeviCell EOS System with 1 year warranty
1000025	LeviCell EOS-4 TEC Core with 1 year warranty
1000026	LeviCell EOS-4 TEC Core with 1 year warranty
<b>Consumables</b>	
1002010	LeviCell Cartridge S2.3 (10 pack)
1002012	LeviCell Cartridge S2.3, sterile (10 pack)
1002101	LeviCell EOS-4 Cartridge (10 pack)
1002102	LeviCell EOS-4 Cartridge, sterile (10 pack)
<b>Reagents</b>	
1003001	Levitation Agent (10 Reaction Kit)
1003002	Levitation Agent (40 Reaction Kit)

Part Numbers	Product Description
<b>LeviSelect Kits</b>	
1004001	LeviSelect Human CD45 Depletion Kit (10 rxn)
1004003	LeviSelect Mouse CD45 Depletion Kit (10 rxn)
1004005	LeviSelect Human CD3 T Cell Enrichment Kit (10 rxn)
1004007	LeviSelect Human CD4 T Cell Enrichment Kit (10 rxn)
1004009	LeviSelect Human CD8 T Cell Enrichment Kit (10 rxn)
1004011	LeviSelect Human B Cell Enrichment Kit (10 rxn)
1004019	LeviSelect Mouse B Cell Enrichment Kit A (10 rxn)
<b>LeviPrep Kits</b>	
1005001	LeviPrep™ Mouse Tissue Dissociation Kit

## Contact Us

sales@levitasbio.com  
levitasbio.com  
650.204.1185



For research use only. Not for use in diagnostic procedures. © 2023 LevitasBio, Inc. All rights reserved. LeviCell, LeviSelect, LeviPrep, LeviMetrics and the LevitasBio logo are trademarks of LevitasBio, Inc. sales@levitasbio.com 90-00250A-0823

**LEVITAS** BIO®

Industriestrasse 12  
CH-6210 Sursee  
mail@witec.ch  
T 041 250 53 57

**witec**  
experts in life science products