

Ordering information

Product	Order number
labcycler Gradient Without block	011-101
labcycler Basic Without block	011-103
Inter System Copy Cable	011-702
Gradient Upgrade (Only for labcycler Basic)	011-801
Thermoblock 384 For microtiterplates 384-well	012-101
Thermoblock 48 For reaction tubes of 0.5 ml	012-102
Thermoblock 96 For reaction tubes of 0.2 ml and microtiterplates 96-well	012-103
Triple Block Without passive lid	012-104
Passive Lid 3 lids are necessary for Triple Block application	012-201
Sealing Pad for Thermoblock 384	012-701

Technical data at a glance

	Device:	labcycler Basic & labcycler Gradient
	Line voltage:	85 V to 265 V without switching, 50 to 60 Hz
	Power Consumption:	Maximum 350 W, standby 25 W
	Loudness:	ldle 38 dBA, typical run 44 dBA, maximum run 48 dBA
	Interfaces:	RS232
	Heated lid:	Electrically moving, temperature and pressure programmable
	Pressure:	Programmable from 30 to 120 Newton
	Dimension:	Length = 444 mm Width = 251 mm Height: lid closed = 201 mm, lid open = 347 mm
	Weight:	11.5 kg
	Display:	TFT illuminated colour display ¼ VGA, 5.7" diagonal, 320 x 240 = 76800 pixel touchscreen
	Keyboard:	Numeric silicone keys Virtual keys on the touch screen depending on the context
	Languages:	English, German
	Programs:	680 5-step-programs, or at least 3000 steps The last 16 program runs can be displayed any time.
	Password Protection:	Individual for groups, persons, folders and programs

Blocks:	Thermoblock 48, 96, 384 and Triple Block
Temperature:	- 5.0 °C to 99.9 °C
Uniformity:	± 0.25 °C at 55 °C, ± 0.40 °C at 95 °C
Control accuracy:	± 0.01 °C
Ramp rate:	0.001 °C/s to 5.0 °C/s
De(In)crements:	Temperature \pm 9.99 °C Time \pm 99.99 seconds
Format:	Thermoblock 48 (48-wells, 0.5 ml single tubes) Thermoblock 96 (96-wells, 0.2 ml single tubes, stripes & microtiterplates) and Thermoblock 384 (384-wells, microtiterplates), electroformed gold plated silver, gradient capable (40 °C, ± 20 °C between the narrow sides of the block) heating rate: 4.2 °C/s, cooling rate: 3.6 °C/s
Format:	Triple Block , 3 x 21 wells, anodised aluminium, 3 passive lids, separately controllable, 0.2 ml single tubes, not gradient capable, heating rate: 2.5 °C/s, cooling rate: 2.2 °C/s 3 different PCR processes at the same time

in one device!

Distributor network worldwide.

SensoQuest Biomedizinische Elektronik GmbH Hannah-Vogt-Straße 1 · 37085 Göttingen · Germany

SENS@UEST Biomedical Electronics

Tel. Sales: +49 551 2503244 · +49 176 66646603 Tel. Technique: +49 551 389195-23 · Fax: -24

E-Mail: info@sensoquest.de · www.sensoquest.com

SensoQuest develops and produces thermocyclers which are sold by international distributors since 2005. The team of physicists, engineers, and biologists is very successful with 20 years of experience in the biomedical market. The company currently has the smallest and most versatile Triple Block system worldwide, as well as the only 384-well silver block.

T 041 250 53 57





labcycler

The SensoQuest team has been developing and making thermocyclers since 1990. After all we thought it was time for a new generation, which we came out with in 2005.

The labcycler features a truly intuitive user interface with a coloured touchscreen, a nice design and solid construction. All that comes with a unique block changing system, giving full flexibility for present and future applications. A choice of three **gold plated silver blocks** was designed for high speed, yet low energy consumption and good temperature uniformity. These are complemented by the Triple Block, which lets you run three independent processes on one machine.

Sustainability and good value were prime considerations. The peltier elements were tested to **600,000 cycles without any failures**, giving at least 20 years of lifetime even under the harshest conditions. The silver blocks are electroformed for lowest heat capacity and best heat conductance. This allows high speed with a maximum power of only 350 Watts. The average during a typical run is less than 150 Watts. The result is good performance with **low energy consumption**, low carbon dioxide footprint, less heat in the lab and, last but not least, less noise from the cooling fans.

Precision is further enhanced by a **6-zone temperature regulation** that corrects for any differences between the 6 peltier elements. Each block has its own processor with a continuously self-calibrating temperature measuring circuitry. Indefinite cooling at 4 °C is of course possible, the blocks even go down to -5 °C.

Although the user interface is quite self-explanatory, a context sensitive online help function further assists you, making the manual a rarely used item.

Programs can be copied between two labcyclers via a cable, making it easy to keep several of them "in line".

Of course there is an **automatic restart** after a failure of the power line. The program will continue with the last denaturation step to prevent false annealing.



Triple Block 3 x 21


- Material: anodised aluminium
- · Thermal conductivity: 237 W/mK
- · Heating rate: 2.5 °C/s · Cooling rate: 2.2 °C/s
- 3 independent PCR-runs
- 3 x 21 wells for 0.2 ml caps
- Minimum volume of reaction: 10 µl
- Protection against condensation by 3 Passive Lids
- Separate and parallel monitoring of all blocks







The labcycler has a **TFT display with a touchscreen** featuring alphanumeric and function keys. Familiar buttons and icons enable an intuitive use. The interface "speaks" English and German.

Graphic monitoring allows tracking of the PCR process for single and Triple Blocks. The Triple Block system is displayed with the TFT touchscreen separated in three parts.







Thermoblocks

With the unique **quick block changing system**, a block change takes one hand and ten seconds.

All thermoblocks have their own processor with **6 separately controlled peltier elements** for extraordinary temperature uniformity at high heating and cooling rates.

The temperature measuring system is entirely in the block and continuously **self-calibrating**, ensuring precise and identical operation of a block in any machine.

Thermoblock 48 Thermoblock 384 Thermoblock **96** Material: Electroformed gold plated silver Thermal conductivity: 429 W/mK Heating rate 4.2 °C/s · Cooling rate 3.6 °C/s 48 well block 96 well block 384 well block 24 zone gradient 8 zone gradient 12 zone gradient 0.5 ml tubes 0.2 ml tubes Gradient capable: 40 °C, ± 20 °C from the left to the right 96 Well microtiterplates 384 Well microtiterplates Minimum reaction volumina 20 µl 10 µl 3 µl

Automatic Lid

The heated lid is controlled by an electric motor. Pressure and temperature are fully programmable.

It quickly reaches its uniform temperature through high power.

During a programmed or manual pause the lid comes up to give access to the probes for **hotstart-procedures**. The temperature and force of the lid can be preselected for each program.

