

Pipette Tips Compatibility Chart

Note: This chart is both suitable for Filter Pipette Tips and non-filtered Tips

Brand	Series	Channel	Volume(μL)	10μL	10μL Extended	20μL	100μL	200μL	300μL	1000μL	1250μL
Axygen	Axypet	Single Channel	0.1-2	✓	✓						
			0.5-10	✓	✓						
			2-20			✓	✓	✓	✓		
			5-50			✓	✓	✓	✓		
			10-100			✓	✓	✓	✓		
			20-200			✓	✓	✓	✓		
			100-1000							✓	✓
		Multi-Channel	1-10	✓	✓						
			5-50			✓	✓	✓	✓		
			20-200			✓	✓	✓	✓		
Biohit	Proline	Single Channel	0.1-2.5	✓	✓						
			0.2-10	✓	✓						
			0.5-10	✓	✓						
			5-50			✓	✓	✓	✓		
			5-100			✓	✓	✓	✓		
			10-250			✓	✓	✓	✓		
			10-500							✓	✓
			20-200			✓	✓	✓	✓		
			50-1000							✓	✓
			100-1000							✓	✓
		Multi-Channel	0.2-10	✓	✓						
			5-100			✓	✓	✓	✓		
			25-250			✓	✓	✓	✓		
			50-250			✓	✓	✓	✓		
			5-100			✓	✓	✓	✓		
	e-Line	Single Channel	0.1-10	✓	✓						
			0.2-10	✓	✓						
			5-120			✓	✓	✓	✓		
			10-300			✓	✓	✓	✓		
			50-1000							✓	✓
		Multi-Channel	0.2-10	✓	✓						
			5-120			✓	✓	✓	✓		
	m-Line	Single Channel	10-300			✓	✓	✓	✓		
			0.1-3	✓	✓						
			0.5-10	✓	✓						
			10-100			✓	✓	✓	✓		
		Multi-Channel	20-100			✓	✓	✓	✓		
			100-1000							✓	✓
			0.5-10	✓	✓						
			5-100			✓	✓	✓	✓		
BrandTech	Transferpett e	Single Channel	30-300			✓	✓	✓	✓		
			0.1-1	✓	✓						
			0.5-10	✓	✓						
			2-20			✓	✓	✓	✓		
			5-50			✓	✓	✓	✓		
			10-100			✓	✓	✓	✓		
			20-200			✓	✓	✓	✓		
			25-250			✓	✓	✓	✓		
		Multi-Channel	100-1000							✓	✓
			0.5-10	✓	✓						
			2-20			✓	✓	✓	✓		
			2.5-25			✓	✓	✓	✓		
			5-50			✓	✓	✓	✓		
			10-100			✓	✓	✓	✓		
			20-200			✓	✓	✓	✓		
			30-300			✓	✓	✓	✓		

Pipette Tips Compatibility Chart

Note: This chart is both suitable for Filter Pipette Tips and non-filtered Tips

Brand	Series	Channel	Volume(μL)	10μL	10μL Extended	20μL	100μL	200μL	300μL	1000μL	1250μL
Capp Aero		Single /Multi-Channel	0.2-2	✓	✓						
			0.5-10	✓	✓						
			5-50			✓	✓	✓	✓		
			10-100			✓	✓	✓	✓		
			25-200			✓	✓	✓	✓		
			50-300			✓	✓	✓	✓		
			100-1000							✓	✓
Eppendorf	Eletronic	Single /Multi-Channel	0.5-10	✓	✓						
			5-100			✓	✓	✓	✓		
			2-200			✓	✓	✓	✓		
			50-1000							✓	✓
	Reference	Single /Multi-Channel	0.1-2.5	✓	✓						
			0.1-10	✓	✓						
			0.5-10	✓	✓						
			2-20			✓	✓	✓	✓		
			10-100			✓	✓	✓	✓		
			20-200			✓	✓	✓	✓		
			50-200			✓	✓	✓	✓		
			100-1000							✓	✓
	Research	Single Channel	0.1-2.5	✓	✓						
			0.5-10	✓	✓						
			2-20			✓	✓	✓	✓		
			10-100			✓	✓	✓	✓		
			20-200			✓	✓	✓	✓		
			100-1000							✓	✓
		Multi-Channel	0.5-10	✓	✓						
			5-100			✓	✓	✓	✓		
			10-100			✓	✓	✓	✓		
			50-100			✓	✓	✓	✓		
			20-300			✓	✓	✓	✓		
			30-300			✓	✓	✓	✓		
			50-1200							✓	✓
			20-300			✓	✓	✓	✓		
	Research Pro	Single Channel	50-1000							✓	✓
		Multi-Channel	0.5-10	✓	✓						
			20-300	✓	✓						
	Research Plus	Single /Multi-Channel	100-1000							✓	✓
Gilson Pipetman	Pipetman Classic	Single Channel	P2,0,1-2	✓	✓						
			P2,0,2-2	✓	✓						
			P10,0,2-10	✓	✓						
			P20,2-20			✓	✓	✓	✓		
			P100,10-100			✓	✓	✓	✓		
			P200,20-200			✓	✓	✓	✓		
			P1000,100-1000							✓	✓
	Pipetman Neo	Single Channel	P2N,0,2-2	✓	✓						
			P10N,1-10	✓	✓						
			P20N,2-20			✓	✓	✓	✓		
			P100N,10-100			✓	✓	✓	✓		
			P100N,20-100			✓	✓	✓	✓		
			P200N,20-200			✓	✓	✓	✓		
			P8x20N,2-20			✓	✓	✓	✓		
	Pipetman Ultra	Single Channel	U200,20-200			✓	✓	✓	✓		
		Multi-Channel	P200-M8,20-200			✓	✓	✓	✓		

Pipette Tips Compatibility Chart

Note: This chart is both suitable for Filter Pipette Tips and non-filtered Tips

Brand	Series	Channel	Volume(μL)	10μL	10μL Extended	20μL	100μL	200μL	300μL	1000μL	1250μL
Labnet	Labpette Discovery or Biopette™ Plus	Single Channel	0.1-2	✓	✓						
			0.5-10	✓	✓						
			2-20			✓	✓	✓	✓		
			10-100			✓	✓	✓	✓		
			20-200			✓	✓	✓	✓		
			100-1000							✓	✓
		Multi-Channel	1-10	✓	✓						
			5-50			✓	✓	✓	✓		
			20-200			✓	✓	✓	✓		
			50-300			✓	✓	✓	✓		
	BioPette E	Single Channel	0.5-10	✓	✓						
			2-20	✓	✓						
			10-200			✓	✓	✓	✓		
			100-1000							✓	✓
		Multi-Channel	0.5-10	✓	✓						
			2-20	✓	✓						
			10-200			✓	✓	✓	✓		
Matrix	Impact/ Impact ²	Single /Multi-Channel	2-125			✓	✓	✓	✓		
Nichryo/ Oxford	Benchmate	Single Channel	0.5-10	✓	✓						
			2-20			✓	✓	✓	✓		
			10-50			✓	✓	✓	✓		
Rainin	EDP3	Single Channel	E-MIC-10,0.5-10	✓	✓						
			E2-25,2.5-25			✓	✓	✓	✓		
			E2-100,10-100			✓	✓	✓	✓		
			E2-250,25-250			✓	✓	✓	✓		
			E2-1000,100-1000							✓	✓
		Multi-Channel	EP-M8-10,0.5-10	✓	✓						
			EP-M8-250,10-250			✓	✓	✓	✓		
	Latch Mode	Single Channel	R2,0.1-2	✓	✓						
			R10,0.2-10	✓	✓						
			R20,2-20			✓	✓	✓	✓		
			R100,10-100			✓	✓	✓	✓		
			R200,20-200			✓	✓	✓	✓		
			R1000,100-1000							✓	✓
Socorex	Callbra/ Accura	Single Channel	0.1-2	✓	✓						
			0.2-2	✓	✓						
			0.5-10	✓	✓						
			1-10	✓	✓						
			2-20			✓	✓	✓	✓		
			10-100			✓	✓	✓	✓		
			20-200			✓	✓	✓	✓		
			100-1000							✓	✓
		Multi-Channel	1-10	✓	✓						
			5-50			✓	✓	✓	✓		
			10-100			✓	✓	✓	✓		
			20-200			✓	✓	✓	✓		
			50-200			✓	✓	✓	✓		

Pipette Tips Compatibility Chart

Note: This chart is both suitable for Filter Pipette Tips and non-filtered Tips

Brand	Series	Channel	Volume(μL)	10μL	10μL Extended	20μL	100μL	200μL	300μL	1000μL	1250μL
Thermo Scientific	Finnpipette	Single Channel	0.2-2	✓	✓						
			0.3-3	✓	✓						
			0.2-10	✓	✓						
			0.5-5	✓	✓						
			0.5-10	✓	✓						
			1-10	✓	✓						
			5-40			✓	✓	✓	✓		
			5-50			✓	✓	✓	✓		
			10-100			✓	✓	✓	✓		
			40-200			✓	✓	✓	✓		
			30-300			✓	✓	✓	✓		
			100-1000							✓	✓
			200-1000							✓	✓
		Multi-Channel	1-10	✓	✓						
			5-50			✓	✓	✓	✓		
			50-300			✓	✓	✓	✓		
	Finnpipette Digital	Single Channel	40-200			✓	✓	✓	✓		
	Finnpipette electronic	Multi-Channel	1-10	✓	✓						
			30-300			✓	✓	✓	✓		
			100-1000							✓	✓
VWR	VWR and VWR Signature	Single Channel	0.1-2	✓	✓						
			0.5-10	✓	✓						
			2-20			✓	✓	✓	✓		
			5-50			✓	✓	✓	✓		
			10-100			✓	✓	✓	✓		
			20-200			✓	✓	✓	✓		
			50-250			✓	✓	✓	✓		
			100-1000							✓	✓
		Multi-Channel	1-10	✓	✓						
			5-50			✓	✓	✓	✓		
			20-200			✓	✓	✓	✓		
			50-300			✓	✓	✓	✓		