Incubation Trays

Convenient, easy to use trays for staining and washing gels and membranes

These lidded **incubation trays** have a perfectly smooth interior surface, making them ideal for staining and washing electrophoresis gels and membranes. Trays are manufactured to prevent adsorption of antibodies. Three sizes are available to hold different volumes or gel sizes, so you can find the right tray for your application.



Advantages

- Smooth interior protects membranes and gels from scratches
- Lidded to protect experiments from dust or debris
- Non-binding surfaces no antibody loss when used with common immunoassay buffers
- Convenient multiple sizes available for different applications

Ordering Information		
Catalog Number	Product	Size
L-07031-005	Incubation Tray, red, 5x7 cm	5/pack
L-07032-005	Incubation Tray, black, 6x9 cm	5/pack
L-07033-005	Incubation Tray, blue, 9x11 cm	5/pack
Related Products		
K-12045-C20	WesternBright™ ECL Western Blotting HRP Substrate Trial kit size	20 ml
K-12045-D20	WesternBright™ ECL Western Blotting HRP Substrate (for 2000 cm² membrane)	200 ml
K-12045-D50	WesternBright™ ECL Western Blotting HRP Substrate (for 5000 cm² membrane)	500 ml
R-05071-500	Goat-anti-mouse HRP-conjugated secondary antibody	500 µl
R-05072-500	Goat-anti-rabbit HRP-conjugated secondary antibody	500 µl
L-08001-010	Low-Fluorescence PVDF Transfer Membrane 7x9 cm	10 sheets
L-08002-010	Nitrocellulose Transfer Membrane 0.45 µm 8x10 cm	10 sheets
L-08003-010	Nitrocellulose Transfer Membrane 0.22 µm 8x10 cm	10 sheets

Copyright © 2010 Advansta. All rights reserved. The Advansta logo is a registered trademark of the Company. WesternBright™ is a trademark of the Company. All other trademarks, service marks and tradenames appearing in this brochure are the property of their respective owners.

Advansta Corporation

1455 Adams Drive, Ste. 1160 | Menlo Park, CA 94025 Tel: 650.325.1980 | Fax: 650.325.1904 | Email: sales@advansta.com Product information: www.advansta.com/Incubation_trays

