




C-Slide

— *New Name
of countless slide*




For used with Automated cell counter



C-Slide

Cell counting chamber slides

- Do not need coverslips
- Eliminate the need of washing and reuse
- Reduce exposures to the infectious samples and contamination
- Guaranteed reproducibility and reliability
- Made of sturdy and strong quartz grade optical plastic



NanoEnTek

The developer & manufacturer
of Countess and Countess slide.



Specifications

Sample volume : 10 μ L

Staining method : Trypan blue stain

Dimension : 75 (L) x 25 (W) x 1.8 (H) mm

For used with Automated cell counter



Countess



Countess II

Ordering information

Catalog Number	Product Name
CS-050	C-Slide, Cell counting chamber slides 50 slides (100 tests)

NanoEnTek

e-mail sales@nanoentek.com
website www.nanoentek.com

NanoEnTek, Inc.

851-14, Seohae-ro, Paltan-myeon, Hwaseong-si,
Gyeonggi-do, 18531, Korea
Tel : +82-2-6220-7940 / Fax : +82-2-6220-7999

Industriestrasse 12
CH-6210 Sursee
mail@witec.ch
T 041 250 53 57

witec
experts in life science products