

## Labnet Mini Centrifuges and Slide Spinner

The Labnet Mini Microcentrifuges operate at a fixed speed of 6,000 rpm (2,000 x g) and were designed for quick spin downs of micro-samples. Operation is made simpler and more convenient with the design of a quick-release rotor system and electronic brake. Once samples are loaded and the lid is closed, the rotor rapidly accelerates to 6,000 rpm. This speed range is ideal for bringing small droplets to the bottom of the tubes for micro-filtrations or as basic separations. Pressing the lid release button will activate the electronic brake, which brings the rotor to a fast stop. The mini microcentrifuge can accommodate up to 8 microcentrifuge tubes or 4 PCR strip tubes. Smaller samples such as 0.2, 0.25, and 0.5 mL can also be processed by using the included tube adapters.

The Labnet Slide Spinner is designed for drying microarray slides before scanning. It dries two slides in as little as 10 seconds for quality imaging and screening results. The slides are placed into cassettes that catch all liquid without splashing. Cassettes are reusable.



Mini  
Microcentrifuge



Slide Spinner

SPECIFICATIONS	Mini Centrifuge	Slide Spinner
Maximum speed/RCF	6,000 rpm/2,000 x g <sup>1</sup>	4,000 rpm/500-850 x g
Maximum capacity	4 x 0.2 mL strips or 8 x 1.5/2.0 mL	2 standard slides
Dimensions (W x D x H)	6.8 x 5.8 x 4.8 in/17.2 x 14.7 x 12.2 cm	5.9 x 5.9 x 4.6 in/15 x 15 x 11.7 cm
Weight	2.65 lb/1.2 kg	1 lb/0.45 kg
Electrical	AC 100 - 240V, 50/60Hz	120V~, 60 Hz or 230V~, 50 Hz

<sup>1</sup>Average max rcf in the strip tube rotor is 1,180 - 2,000 x g

CAT NO.	DESCRIPTION
C1601-G	Labnet Mini Microcentrifuge, gray lid, 100-240V
C1601-B	Labnet Mini Microcentrifuge, blue lid, 100-240V
C1601-P	Labnet Mini Microcentrifuge, purple lid, 100-240V
C1601-R	Labnet Mini Microcentrifuge, red lid, 100-240V
6770-RT	Replacement rotor, 8 place tube rotor
6770-RTS	Replacement rotor, 4 place PCR strip rotor
6770-RTA.5	Individual adapters for 0.5 mL tubes, pk of 8
6770-RTA.2	Individual adapters for 0.2 mL tubes, pk of 8
6770-RTA.25	Individual adapters for 0.25 mL tubes, pk of 8
401025	Power Adapter, 24V
C1303-T*	Slide Spinner with two slide drying cassettes, 120V
C1303-SC	Extra cassettes for Slide Spinner centrifuge, pk of 2