

Performance Parameter

Pipette	Pipette Type	Fixed 8 channel	
	Pipette Range	1 μ L~200 μ L	
	Pipette CV	1 μ L : <5%	200 μ L : <1%
	Pipette Accuracy	1 μ L : < \pm 10%	200 μ L : < \pm 1%
PCR	PCR Temp. Range	4~99 $^{\circ}$ C	
	PCR Temp. Accuracy	\pm 0.3 $^{\circ}$ C at 55 $^{\circ}$ C	
	PCR Temp. Uniformity	\pm 0.2 $^{\circ}$ C at 72 $^{\circ}$ C	
Temperature Control Module	Temperature Control Module Temp. Range	4~90 $^{\circ}$ C	
	Temperature Control Module Temp. Accuracy	\pm 1 $^{\circ}$ C at 55 $^{\circ}$ C	
	Temperature Control Module Temp. Uniformity	\pm 1 $^{\circ}$ C at 72 $^{\circ}$ C	
Robot Arm Positioning Accuracy	\pm 0.1mm		
HEPA System	99.99%@0.3 μ m filtration efficiency		
UV Lamp	Radiation dose higher than 100 000 μ W.s/cm ²		
Operating Environment Requirements	Temperature	19 $^{\circ}$ C~25 $^{\circ}$ C	
	Relative Humidity	20%RH~80%RH, non-condensing	
	Atmospheric Pressure	80~106kPa	
	Waterproof Rating	IPX0	
Power Type	100~240V, 50/60HZ, 1600VA		
Dimensions	Length 780mm,Width 725mm,Height 777mm		
Net Weight	130kg		



Design for Massively Parallel Sequencing Library Preparation

Automated

Integrated

Stabilized

Customized



DNA Sequencing Library Preparation System

MGISP-100



Ordering Information

Model	Product Name	Intended Market	Item NO.
MGISP-100RS	DNA Sequencing Library Preparation System MGISP-100RS	CE RUO	900-000206-00
MGISP-100	DNA Sequencing Library Preparation System MGISP-100	CE IVD	900-000207-00

MGI Tech Co., Ltd.

Building 11, Beishan Industrial Zone, Yantian District, Shenzhen

MGI-service@mgi-tech.com en.mgi-tech.com +86 4000-688-114

Copyright Disclaimer

The copyright of this brochure is solely owned by MGI Tech Co. Ltd.. The information included in this brochure or part of, including but not limited to interior design, cover design and icons, is strictly forbidden to be reproduced or transmitted in any form, by any means (e.g. electronic, photocopying, recording, translating or otherwise) without the prior written permission by MGI Tech Co., Ltd.. All the trademarks or icons in the brochure are the intellectual property of MGI Tech Co., Ltd. and their respective producers.

Version: February 2024 | MGPC2002202

Information in this brochure is updated to [02/23/2024] and only for your reference. In no event shall the brochure be regarded as warranty or commitment made by MGI Tech Co., Ltd. All rights and obligations shall subject to final executed agreement.



MGI Tech Co., Ltd.


MGISP-100

MGISP-100 Automated Sample Preparation System is an automated liquid handling workstation specialized for MPS library preparation, it can process samples in batches, eliminate operators from repetitive processing procedures, increase the stability of library preparation, reduce the total costs, and significantly enhance laboratory overall efficiency.

8 channel Pipette

Accurate positioning, fast movement, high efficiency

Pipette Type	Fixed 8 channel pipette
Pipetting Range	1 μ L-200 μ L
Pipetting Precision	1 μ L : <5% 200 μ L: <1%
Pipetting Accuracy	1 μ L: < \pm 10% 200 μ L: < \pm 1%
Positioning Precision of Robotic Arm	\pm 0.1mm



HEPA system

Positive pressure
99.99%@0.3 μ m filtration efficiency

UV Disinfection

Integrated 2 UV lamps, ensure that each space of the instrument can obtain the same intensity of exposure

Radiation dose
higher than 100 000 μ W.s/cm²

Automated PCR Machine

Precise temperature control

Temperature range	4~99 $^{\circ}$ C
Temperature accuracy	\pm 0.3 $^{\circ}$ C at 55 $^{\circ}$ C
Temperature uniformity	\pm 0.2 $^{\circ}$ C at 72 $^{\circ}$ C

Automatic up/down Magnetic Rack

The magnetic rack can be up/down automatically. Avoid the mechanical gripper. Designed for more stable operation

Temperature Control Module


Accurate temperature control, range from 4 $^{\circ}$ C-90 $^{\circ}$ C

Temperature range	4~90 $^{\circ}$ C
Temperature accuracy	\pm 1 $^{\circ}$ C at 55 $^{\circ}$ C
Temperature uniformity	\pm 1 $^{\circ}$ C at 72 $^{\circ}$ C

Product Feature

- Automated** Easy-to-use, automated 90% steps, reduce repetitive processing procedures.
- Integrated** Compact design,integrated liquid handling, automated magnetic frame and PCR machine into less than 0.5m³.
- Stabilized** Integrated UV lamp and HEPA system and minimum 1uL precise pipetting supported.
- Customized** Support customized development according to different needs and applications, graphical drag-and-drop script editor, which is compatible with Python diverse script development.

WD graphical drag-and-drop script editor



Automated Package for MPS Library Preparation

